



NISUS®

NEGATIVE PRESSURE WOUND THERAPY PUMP



Latch design provides streamlined, one-handed canister release and removal



User friendly interface with easy to read LCD screen

The Cork Nisus Negative Pressure Wound Therapy Pump is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.

NISUS NPWT PUMP FEATURES

The Cork Medical Products Nisus Negative Pressure Wound Therapy Pump is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- Recommended for patients with chronic, acute, traumatic, sub-acute and dehiscent wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, and venous ulcers, flaps, and grafts
- FDA approved for use in hospitals, long term care facilities, and home-care settings
- Light-weight and portable, weighing less than 2 lbs. with up to 24 hours battery life (can extend to 48 hours with sleep mode)
- Sleep feature technology allows user to dim light for longer lasting battery life and discreetness
- Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage
- Multiple features for clinical enhancement such as: patient runtime, alarm log, lock restriction, Variable Intermittent and Continuous Therapy Modes

COMPONENT SPECS

The Nisus Negative Pressure Wound Therapy Pump is part of a complete wound care system. Two essential components that actively work together to promote wound healing are the Nisus NPWT Canisters and the Cork NPWT Wound Dressing Kit.

Please use under the advisement of a physician. Carefully read all therapy and device instructions and safety information prior to use. Use only Cork Medical Products components with the Nisus Negative Pressure Wound Therapy Pump.

- Battery charger and carrying bag included
- Medicare HCPCS Code: E2402

Model No.

CMPP-100
CPC-250
CPC-500
WDK-250
WDK-300
WDK-400
CLEANSE-KIT-A-18
NIS-BRDR-A-16
FOAM-WHITE-A-16
CSP-100-ES-EXT
CWY-100

Component

Nisus NPWT Pump
Nisus NPWT 250mL Canister
Nisus NPWT 500mL Canister
NPWT Wound Dressing Kit (Small)
NPWT Wound Dressing Kit (Medium)
NPWT Wound Dressing Kit (Large)
NPWT Cleanse Therapy Kit
Bridge Foam
HydroLite White Foam
24" Extension Tubing
Y-Connector



AGF ENTERPRISE
Real Solutions For Real People With Real Problems

1-866-216-1186

Revision 2.1 2018_1226